




# Oxford Microplasty Instrumentation

CSAO Conference 46<sup>th</sup> Annual  
September 13<sup>th</sup>-15<sup>th</sup>, 2015



# Oxford Microplasty Instrumentation

Proven, safe and reproducible results




# Topics

- Instrumentation
  - Set Overview
- Surgical Technique
  - Key Steps
    - Highlighting Key Instruments
    - Improvements vs Phase 3 technique




# Instrumentation





# Instrumentation Overview

**Phase 1**  
1976 - 1988



- Problems with balancing due to femoral prep
- A lot of eyeballing

**Phase 2**  
1988 - 1998



- Introduction of milling
- Improving reproducibility

**Phase 3**  
1998 - 2011



- Improving milling technique
- Continued focus on reproducibility

**Microplasty**  
2011 – Present...



- Focus on reproducibility

# Instrumentation: Set Overview

## *Tibia 1*


### **Top Tray**

- Tibial Templates and Trials
- Nails, drill bits, puller



### **Bottom Tray**

- Tibial Impactor, Inserter
- Tibial Groove Cutters
- Slap Hammer
- T-Handle
- Bearing Inserter/Extractor




# Instrumentation: Set Overview

## Tibia 2


### Top Tray

- Tibial Resection
- IM Link
- IM Rod Pusher



### Bottom Tray

- |                 |                       |
|-----------------|-----------------------|
| • IM Rods       | • Drills              |
| • Anterior Mill | • Bone Collar Remover |
| • Impactor      | • 5mm Awl             |
| • Spigots       |                       |
| • Toffee Hammer |                       |




# Instrumentation: Set Overview


## *Size Specific*

### **Specific for Femoral Size**


- Femoral Trials
- Bearing Trials
- Femoral Drill Guide
- Femoral Resection Guide
- Anti-Impingement Guide
- Feeler Gauges
- Mill




# Overview of Size Specific Trays




**XSM**




**SM**



**MD**



**LG**



**XLG**



# Surgical steps

highlights of Oxford Microplasty



# Oxford Microplasty

- Design Goal: Simple, Easy, Accurate, Reproducible
- No change in indications
- Enhancement in tibial resection instruments
- IM linked femoral preparation
- Anti-impingement instrumentation
- Trial bearing inserts by hand

# Indications/Contraindications

## Indications

- Antero-medial OA
- Intact ACL
- Full Thickness Lateral
- Correctable Varus
- Fixed Flexion Deformity (FFD) less than 15°

## Not Contraindication

- Patello-femoral Joint
- Obesity
- Age
- Activity

1 in 4 OA Knees!

# Tibial Preparation

## Goals

- Minimize resection depth
- Maximize surface area
- Decrease re-cutting
  - Offers shims for recutting in +/- 2mm increments
- Appropriate slope



# Tibial Preparation


- Tibial Slope
  - Parallel to tibial crest
  - 7 degrees
- New ankle positioning
  - Helps avoid overslope
- Pinned to tibia
  - Single nail suffices



# Tibial Preparation

- Sizing spoon placed under posterior medial femoral condyle
- 1,2,3 mm spoons (1mm spoon is default)
- Ligament tension achieved
- Clamped to tibial guide with G-Clamp





# Tibial Preparation


- Vertical cut important
- M/L position
  - Adjacent to ACL footprint
- Rotational position
  - ASIS/Flexion plane
- Depth
  - DON'T LIFT HAND



# Tibial Preparation

- “Curly Whirly” inserted
  - Protect the MCL
- Flat cut
- Slotted guide available
- Shim option
  - Additional 2mm





# Femoral Preparation

- IM rod cannulated
  - Flexion/extension
  - Varus / valgus
- Starting position critical
- Yellow pusher
- Mark central third




# Femoral Preparation


- Femoral drill guide set
- Inserted and linked
- Position of link important
- Must be seated
- M/L position confirmed



# Femoral Preparation

- “Flat arm” of link in IM hole
- Curved arm in femoral drill guide
- Foot of guide against tibial resection
  - In contact with femur
- Central third confirmed
- Correct hole, not femoral drill
  - Link out of way of drill hole
  - Places drill holes 10° flexed, and 7° valgus






# Femoral Preparation


- New curved slotted posterior cutting guide
- Suitable for Oxford oscillating sawblade
- Allows blade to be flexed and for saw throw
- Goal: to prevent over/under resection





# Femoral Preparation

- Milling
  - Curve off femur with 0 spigot



# Gap Assessment

- Use single peg femoral trial to balance
- New tapered feeler gauges
- Less soft tissues interference
- Gaps determined
- Standard milling to balance
- 1,2,3 mm plastic feeler gauges




# Anti-Impingement

- Decreases early and late complications
- Replaces traditional “Freehand technique”
- New instruments help guide bone removal



# Anti-Impingement

- Anterior mill
- Size specific
- Similar to a calcar reamer
- Spring loaded plunger
- Extend the knee
- Care to watch tibia...




# Final Trialing

- Insert tibial tray trial and single peg femoral trial
- Trial bearing inserted by hand
- Better “feel” of tension
- Assess tracking
- Check impingement



# Trial Bearing Removal



# Cementing the Components

- Cementing technique follows current Phase 3 Oxford Knee surgical technique
- Two staged
  1. Tibia
    - Small amount on tibia
    - Use tibial impactor
    - Compress at 45 deg with trial femur
  2. Femur
    - Cement into large hole
    - Concave surface of femoral component
    - Impact at 45 deg to the long axis
    - Compress at 45 deg with feeler gauge





# Cementless

Differences vs Cemented

# Keel Preparation

Use cementless tooth brush saw  
through tibial template

**Cemented**

Use cemented tibial groove cutter


**Cementless**

Option to use the cementless tibial  
groove cutter



# Implant

- Use cementless tibial inserter
- Two insertion tabs on cementless tibia
- Bring in at an angle, and lightly tap home





# Conclusion

# Oxford Microplasty

*Proven to deliver more accurate and reproducible results*

- 219 Ph3 v 196 OXMP
- Concluded that OXMP delivered:
  - A more accurate alignment of the femoral implants
  - A more accurate alignment of the tibial implants
  - A more accurate resection, with a greater number of thinner 3mm and 4mm bearings

**ARTICLE IN PRESS**

The Journal of Arthroplasty xxx(2014) xxx-xxx

Contents lists available at ScienceDirect

The Journal of Arthroplasty

journal homepage: [www.arthroplastyjournal.org](http://www.arthroplastyjournal.org)

Volume 29 Number 8 December 2014

ISSN 1093-8143

1093-8143/\$ - see front matter © 2014 Elsevier Inc. All rights reserved.

**Radiographic Comparison of Mobile-Bearing Partial Knee Single-Peg Versus Twin-Peg Design**

Jason M. Hurst, MD <sup>a</sup>, Keith R. Berend, MD <sup>a,b,c</sup>, Joanne B. Adams, BFA <sup>a</sup>, Adolph V. Lombardi Jr., MD, FACS <sup>a,b,c</sup>

<sup>a</sup> John Implant Surgeons, Inc., New Albany, Ohio

<sup>b</sup> Department of Orthopaedics, The Ohio State University Wexner Medical Center, Columbus, Ohio

<sup>c</sup> Mount Carmel Health System, New Albany, Ohio

**ARTICLE INFO**

**Article history:**  
Received 25 July 2014  
Accepted 10 October 2014  
Available online xxxx

**Keywords:**  
partial knee arthroplasty  
mobile-bearing component  
single-peg design  
radiographic analysis  
implant position

**ABSTRACT**

The femoral component and patellofemoral instrumentation of a mobile-bearing, unicondylar experimental knee arthroplasty (UKA) was redesigned with an additional peg for enhanced fixation, 15° of extra femoral surface for contact in deep flexion, more rounded profile, better fit into the medial surface, and updated instrumentation based instrumentation. To assess the benefit of these changes, we compared preoperative radiographs of 219 single-peg mobile-bearing partial knee designs (Fig. 1A) with 196 enhanced twin-peg designs (Fig. 1B) and evaluated them with radiolucent markers. Radiographic analysis showed improved and consistent component positioning with the twin-peg design compared with updated instrumentation compared with the single-peg. The enhanced pegs facilitate improved implant position by using the twin-peg component and updated instrumentation. The clear and carry forward our potential. More robust follow-up is imperative.

© 2014 Elsevier Inc. All rights reserved.

Numerous studies have reported favorable results with a modified mobile-bearing partial knee device (UKA) [1–13]. The most recent modifications to the mobile-bearing Oxford Partial Knee (Biomet, Warsaw, IN) were intended to improve the process of implantation and the reproducibility of implant positioning, and enhance the stability and mechanics of the femoral component, as an attempt to improve the reproducibility of implant position, an intraoperative based femoral preparation guide (Fig. 2A and B) was developed (Fig. 3) and 27 of these femoral surfaces (Fig. 1B) were instrumented (Biomet). This new instrumentation was specifically designed to incorporate the newly modified femoral component. The femoral component was redesigned with an additional peg for more secure fixation and 15° of extra femoral surface for contact in deep flexion. Further design modifications have been made to the twin-peg to enhance kinematics and reduce the risk of component wear. In addition, the enhanced twin-peg (Fig. 1B) was developed to fit into the medial, lateral, and distal femoral surfaces (Fig. 1B). To assess the benefit of these instrumentation and design changes, we compared the preoperative radiographs and implantation related complications of the single-peg design (Fig. 1A) with the enhanced twin-peg design (Fig. 1B) implanted with the updated instrumentation.

**Materials and Methods**

A query of our practice's arthroplasty registry revealed 359 patients (413 knees) who signed an IRB-approved general research consent allowing retrospective review, and underwent medial unicondylar partial knee arthroplasty performed with a mobile-bearing device between August 2008 and December 2011. Use of the enhanced twin-peg and updated instrumentation was commenced in July 2011. One hundred eighty-eight knees (223 knees) received the single-peg design (Fig. 1A) and 171 (190 knees) received the enhanced twin-peg (Table 1). The preoperative diagnosis was associated necrosis in 1 knee (enhanced twin-peg group) and osteoarthritis in all others. The groups were well matched in terms of gender, age, body mass index, preoperative ROM, and Knee Society pain and clinical scores. Forty-seven percent of patients were males ( $n = 167$ ) and 53% were females ( $n = 192$ ). Mean age was 63.5 years, mean weight of 150.0 kg, mean height of 58.4 in, mean BMI was 32.2 kg/m<sup>2</sup>, mean range of 17.5–74.7°/m<sup>2</sup>, and mean ROM was 115.8° (SD 10.2, range 20°–135°). Patients in the single-peg group had slightly higher preoperative Knee Society function scores (51.4 versus 52.3, 0–100 possible).

Radiographs of adequate quality for review were available for 219 single-peg and 186 enhanced twin-peg knees. We evaluate anteroposterior (AP) and lateral radiographs of each knee were evaluated by a single observer (MH) according to criteria described in the manufacturer's surgical technique guide (Oxford Partial Knee Microplasty Instrumentation Technique). Figs. 3 and 4 depict these criteria and measurements, and are as follows: position and size of the femoral component relative to the femur with varus/valgus angle <10° versus to >10° values (Fig. 3 criteria A), flexion/extension angle within 15° flexion or extension (Fig. 4 criteria B), midfemoral placement central (Fig. 3 criteria C),

**Author contributions:** Funding and/or materials support of this study was received from Biomet, Inc. Consulting fees, mobility payments, speakers' honoraria, and research support received from Biomet, Inc.; equity payment received from Biomet, Inc.; research support received from Stryker, Alimed, and Panta. The funders had no role in the design of the study and the collection, analysis, and interpretation of data. The article can be found at <http://dx.doi.org/10.1016/j.arth.2014.10.015>.

**Reprint requests:** Jason M. Hurst, MD, Joint Implant Surgeons, Inc., 7277 Southgate Blvd, New Albany, OH 43110, USA.

**http://dx.doi.org/10.1016/j.arth.2014.10.015**

**0883-5400/\$ - see front matter © 2014 Elsevier Inc. All rights reserved.**

**Please cite this article as:** Hurst JM, et al. Radiographic Comparison of Mobile-Bearing Partial Knee Single-Peg Versus Twin-Peg Design. *J Arthroplasty* (2014). <http://dx.doi.org/10.1016/j.arth.2014.10.015>.

## Support Materials

*Currently available on myBiomet*

- Cementless Microplasty surgical video
- Cemented Microplasty surgical video
- Printed and digital surgical technique for Cementless and Cemented Microplasty
- Cemented Microplasty surgical technique animation
- Phase 3 vs Microplasty instrumentation animation

# Summary

- Oxford Microplasty Instruments enhance
  - Tibial prep, resection depth (femoral size)
  - Femoral prep, component alignment precision
  - Impingement avoidance
- Oxford Microplasty Instruments are more accurate with less outliers than Phase 3 instrumentation



# Disclaimer

*This publication and all content, artwork, photographs, names, logos and marks contained in it are protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer, Biomet or its affiliates, unless otherwise indicated. This publication must not be used, copied or reproduced in whole or in part for any purposes other than marketing by Zimmer, Biomet or its authorised representatives. Use for any other purposes is prohibited.*

*Zimmer and Biomet do not practice medicine and do not recommend any particular orthopaedic implant or surgical technique and is not responsible for use on a specific patient. The surgeon who performs any implant procedure is responsible for determining and utilizing the appropriate techniques for implanting prosthesis in each individual patient.*



ZIMMER BIOMET  
Your progress. Our promise.<sup>™</sup>